

## **INSTRUCTIONS FOR USE**

### **CRYOCHAMBERS**

**Cryochambers are used with Cry-Ac's which is a Hand-Held Cryosurgical Device for the controlled dispensing of Liquid Nitrogen.**

#### **1. Operating Instructions**

Place the open edge of the Chamber on the surface of the tumor to be treated.

Make certain that the entire 360-degree edge is sealing on the surface.

A more uniform freeze will be achieved if the Chamber is held in a vertical position, therefore, we have included the Malleable Extension to facilitate this.

If the vertical position is maintained, the pool of liquid nitrogen will create an extremely deep symmetrical freeze.

#### **2. Cautions**

The vent tube should be pointed in a safe direction away from the operator and the patient. As the freezing progresses, a wet liquid level line will appear on the outside of the Chamber. This line can be maintained with intermittent pressing of the trigger on your CRY-AC. It does not have to be filling the Chamber continually. This will require a pressure of 1 to 2 pounds.

**We emphasize the depth of 3cm to 4cm of destruction may be achieved, therefore significant caution should be taken**

#### **3. Cleaning Instructions**

Perform the manual cleaning procedure as follows:

- a. Prepare a bath of ambient temperature utility (tap) water.
- b. Moisten a lint-free cloth in the utility water bath and wipe the cryochamber for 1 minute to remove as much debris as possible. Use a fresh lint-free cloth as necessary. NOTE: Take care not to get water inside the lumens.
- c. Prepare an Enzol bath using lukewarm utility water, (31-40 degrees C), with a recommendation of 1 oz/gallon.

- d. Immerse each cryochamber in detergent solution. NOTE: Take care not to get water inside the lumens.
- e. After the 1-minute soak time but while immersed, wipe each cryochamber thoroughly with a lint-free cloth for at least 1 minute.
- f. Pay particular attention to hard to clean areas such as rough surfaces and joints. NOTE: Take care not to get water inside the lumens.
- g. Using a fresh lint-free cloth moistened with (deionized) water, wipe each cryochamber thoroughly to remove residual detergent.
- h. Pay careful attention to difficult to clean areas such as joints, and rough surfaces.
- i. Use a lint free cloth to thoroughly dry the cryochambers, including any lumens and rough surfaces. NOTE: Be certain the cryochambers are thoroughly dried before use.

#### **4. Sterilization**

1. Use the following recommended validated sterilization parameters: of sterilization.

- Moist heat sterilization with Gravity cycle is the recommended method of sterilization. Gravity displacement cycle is not recommended.
- Vaporized Hydrogen (VHP), Ethylene oxide (EO), gas plasma and dry heat are not recommended sterilization methods for reusable instruments.
- The recommended parameters demonstrate the minimum validated steam sterilization time and temperature required to achieve a  $1.0 \times 10^{-6}$  sterility assurance level (SAL).
- Sterilize the instruments after placing them in a Stainless-Steel Sterilization Tray and wrapping the tray in a double layer of Bioshield Sterilization Wrap using the envelope technique.
- The validated reprocessing instructions are not applicable to trays that include devices not manufactured or distributed by Brymill.

Cycle Time	Temperature	Exposure Time	Dry Time
Gravity	121°C (250°F)	30	15

- Replace the silicone vent tubing prior to use.

CAUTION: U.S. Federal law restricts this device to sale by or on the order of a physician or veterinarian.

If you have any other questions or comments, please do not hesitate to contact us on (800) 777-2796 (USA) or (0) 1256 841045 (UK) or email us at [brymill@brymill.com](mailto:brymill@brymill.com).

-----  
-  
**EU Authorized Representative –**

Brymill Cryogenic Systems (UK) Ltd, 26 Hayes View Drive, Cheslyn Hay Walsall UK W56 7EX  
email: [sales@brymilluk.com](mailto:sales@brymilluk.com)

