



World Leader in Cryosurgical Equipment Since 1966



#1 Choice of Dermatologists!
85% of dermatologists surveyed use Brymill's CRY-AC®
Largest selection of LN₂ cryosurgical spray tips, probes, & accessories

The Brymill Cry-Ac® Hand Held Liquid Nitrogen Delivery System

Brymill is the world's leading manufacturer of handheld, liquid nitrogen cryosurgery devices used to treat a wide range of common skin lesions.

Our 3 units, available in 500ml, 300ml, and 160ml capacity, provide the ultimate in controlled application of Liquid Nitrogen to the lesion.

The No.1 choice of dermatologists, Brymill offers unsurpassed safety, versatility, and control for precise and effective cryosurgery.



- Fingertip trigger for maximum control and ease of viewing
- Patented Safety Autovent maximizes safety by allowing internal pressure to gradually vent as the cap is unscrewed
- Streamlined relief valve maintains a constant operating pressure to ensure consistent and accurate freezing.
- Polypropylene Cover and Collar insulates the user's hand.
- Stainless steel and brass construction for long life
- Sturdy Delrin base for extra stability

CRY-AC® CRY-AC-3® CRY-BABY®

*A set of 4 Apertures (A, B, C & D), a 1" x 20g Straight Spray and a 3" x 20g Bent Spray are included with each CRY-AC® CRY-AC®-3 and CRY-BABY® for more accurate, controlled treatment of lesions of different sizes and types.

(The "C" Aperture will be attached to the Unit)

Model	Description	Height	Weight		Static Holding Time
			Full	Empty	
B-700	CRY-AC® (500ml) 16oz.	11.25in (286mm)	31.6oz (896g)	19.0oz (539g)	10-12 hours
B-800	CRY-AC®-3 (300ml) 10oz.	8.5in (216mm)	25.2oz (714g)	17.3oz (490g)	6-8 hours
B-400	CRY-BABY®(160ml) 5.4oz.	7.75in (197mm)	23.2oz (658g)	16oz (454g)	3-5 hours

Starter Packs

Family Package

MODEL	DESCRIPTION
BRY-1000	Cry-Ac®-3 – 10oz., Withdrawal Tube, 20SC Dewar, Cryoplate (see page 5)
BRY-1001	Cry-Ac® – 16oz., Withdrawal Tube, 20SC Dewar, Cryoplate (see page 5)

Dermatology Package

MODEL	DESCRIPTION
BRY-1002	Cry-Ac®-3 – 10oz., Withdrawal Device, 20 liter Dewar, Cryoplate (see page 5)
BRY-1003	Cry-Ac® – 16oz., Withdrawal Device, 20 liter Dewar, Cryoplate (see page 5)

Additional Packages *NEW*

MODEL	DESCRIPTION
BRY-1004	Cry-Ac®-3 – 10oz., Withdrawal Device, 30 liter Dewar, Cryoplate (see page 5)
BRY-1005	Cry-Ac® – 16oz., Withdrawal Device, 30 liter Dewar, Cryoplate (see page 5)
BRY-1006	Cry-Ac®-3 – 10oz., Withdrawal Tube, 30 liter Dewar, Cryoplate (see page 5)
BRY-1007	Cry-Ac® – 16oz., Withdrawal Tube, 30 liter Dewar, Cryoplate (see page 5)

Liquid Nitrogen Storage & Withdrawal Equipment



Size	Model	Static Holding Time	Weight Empty	Weight Full	Dimensions	
					Overall Height	Outside Diameter
5 Liter	501-5	4 – 5 wks.	8 lbs.	17 lbs.	18.2in. 462mm	8.8in. 222mm
10 Liter	501-10	6 – 8 wks.	13 lbs.	30 lbs.	21.5in. 546mm	10.3in. 260mm
20 Liter	501-20	8 – 12 wks.	19 lbs.	53 lbs.	24.7in. 627mm	14.5in. 368mm
20 Liter	501-20SC	Long Last 200 days	19 lbs.	53 lbs.	24.7in. 627mm	14.5in. 368mm
30 Liter	501-30	14 – 16 wks.	26 lbs.	77 lbs.	24.1in. 611mm	17in. 432mm
50 Liter	501-50	14 – 17 wks.	34 lbs.	123 lbs.	30.7in. 79mm	17in. 432mm



Withdrawal Tube

Model

Description

- 502** Withdrawal Device (*) - For Dewars of 20 liters and over
- 503** Withdrawal Tube (*) - For Dewars of 5 up to 20 liters
- 502-F** Replacement Filter for Withdrawal Device
- 503-F** Replacement Filter for Withdrawal Tube



Withdrawal Device



Roller Base (*)

Recommended for Dewars of 20 liters and over
501-RB-S for 20L
501-RB-L for 30L & 50L



Dewar Cap (*)

Replacement
501-DC-2.15 for 5,10,20L
501-DC-2.5 for 30,50L
501-DC-2 for 20SC older than 2015

Note: When ordering the asterisked items, please specify Dewar manufacturer, model # and capacity



Model

Description

- 605-S** Cryo Gloves – Small
- 605-M** Cryo Gloves – Medium
- 605-L** Cryo Gloves – Large

Mid-arm length: 14 - 15in
 Recommended to protect people withdrawing or dispensing liquid nitrogen from a Storage Dewar

Open Sprays

		
101-AA	102-C	301-E
Model	Description	
101-AA	Advanced Acne Aperture	
102-A	Aperture A (0.04in.)	
102-B	Aperture B (0.031in.)	
102-C	Aperture C (0.022in.)	
102-D	Aperture D (0.016in.)	
102-6	Tip Package - Set of 4 Apertures (A-D), 1" x 20g Straight Spray and 3" x 20g Bent Spray.	
301-E	Super E Tip (0.013in) with back vent	

106-20	
107-16	
107-18	
107-20	
Model	Description
106	Original Straight Spray (not pictured)
106-20	3in. x 20g Straight Spray
107-16	1.0in. x 16g Straight Spray Needle
107-18	1.0in. x 18g Straight Spray Needle
107-20	1.0in. x 20g Straight Spray Needle

Spray Extensions & Adapters

	Model	Description
	103	Original Bent Spray
	103-20	3in. x 20g Bent Spray
	103-20S	1.5in. x 20g Bent Spray

			
301	308	309	310-45 & 310
Model	Description		
301	Back Vent Adapter		
308	Luer Lock Adapter		
309	Malleable Extension		
310	Right Angle Adapter (Recommended for use with Cryo Chambers)		

Cryochambers

	Model	Description
	104-6	6mm Cryochamber
	104-10	10mm Cryochamber
	104-15	15mm Cryochamber
	104-18	18mm Cryochamber

Cryoplate

	Model	Description
	303	Cryoplate
		A transparent lexan plate with 4 various sized conical openings (3, 5, 8 and 10mm). Used with apertures, the cryoplate provides even greater localization of freezing and protection of sensitive areas such as eyes.

Closed Probes

Mini Probes

	Model	Description
	214-1	1mm Mini Probe
	214-2	2mm Mini Probe
	214-3	3mm Mini Probe
	214-4	4mm Mini Probe
	214-5	5mm Mini Probe
	214-6	6mm Mini Probe
	214-SP	Sharp Pointed Mini Probe

Conical Probes

	Model	Description
	216	Sharp-pointed Conical Probe
	203-1	1mm Conical Probe
	203-2	2mm Conical Probe
	203-3	3mm Conical Probe
	203-4	4mm Conical Probe
	203-5	5mm Conical Probe
	203-6	6mm Conical Probe

Ball Probes


	Model	Description
	201-6	6mm Ball Probe
	201-8	8mm Ball Probe
	201-1	1cm Ball Probe
	201-2	2cm Ball Probe

Flat Probes

	Model	Description
	205-1	1cm Flat Probe
	205-2	2cm Flat Probe
	205-2.5	2.5cm Flat Probe
	205-3	3cm Flat Probe
	205-8	8mm Flat Probe

Additional Products

	Model	Description
	601	Cryotote
	Durable insulated bag with long handle for convenient hands-free carrying of CRY-AC® or CRY-AC®-3 off-site	

	<i>Cutaneous Cryosurgery: Principles & Clinical Practices: - Third Edition.</i> The leading textbook that covers essential information for starting and maintaining a cryosurgical practice. Describes diagnoses and cryosurgical treatments for many benign, pre-malignant, and malignant lesions. Extensively illustrated with numerous color photos that can aid in diagnosis and treatment of various dermatological conditions. By Drs. Rodney Dawber, Graham Colver, and Arthur Jackson. Available from your Medical Book Supplier ISBN number: - 1084184-552-3
---	---

	Model	Description
	350	Cryo Tweezers
	The Brymill Cryo Tweezers have been designed for the localized treatment of Skin Tags. By using the Brymill Cryo Tweezers you can control freezing easily and accurately. Supplied with a Delrin Cup Holder and 10 Foam Cups. To reorder more Styrofoam Cups simply Google Dart 4J4	

	Model	Description
	360	Histology Package
	The Brymill Histology Package consists of a 10-inch Malleable Spray Extension and a Wide Spray (Acne) Aperture.	

Quality Assurance Guaranteed

All products designed and manufactured by Brymill Cryogenic Systems are done so in accordance with appropriate US FDA, Health Canada and EU Regulatory requirements, underscoring Brymill's commitment to supplying safe, quality products.



Brymill Offices



Brymill Cryogenic Systems

World Headquarters

105 Windermere Avenue
Ellington, CT 06029
Tel: (860) 875 2460
Email: brymill@brymill.com
Web: www.brymill.com

Fax: (860) 872 2371

California Office

Tel: (619) 925 4796
Email: dmbryne@san.rr.com
Web: www.brymill.com



Brymill Cryogenic Systems – Europe

European Authorized Representative

Rosa Maria Sallent Moya

Calle Jaume Balmes, 59

08830 Sant Boi de Llobregat

SPAIN

Phone: +34 616602892

Email: brymill.internationalsales@gmail.com

Brymill's US Dealers and International Distributors

To find your nearest US Dealer, International Distributor or Repair Center visit our web site at http://www.brymill.com/distributors_repair_centres.html

Why Brymill

- World Leader Since 1966
- #1 Choice of 85% of Dermatologist
- Largest Selection of LN2 Equipment, Tips Probes & Accessories
- Superior Quality, Reliability & Durability
- 3-Year Limited warranty against manufacturing defects

What Dermatologists Say...

"Brymill hand-held units are an indispensable part of my practice. Because of the sizes and design, they are comfortable to hold and easy to use. They function excellently, without clogging, and seem to last forever. The large selection of accessories available makes it possible for me to treat many types of benign, pre-malignant and malignant lesions. I would recommend these units to any physician interested in performing cryosurgery."

*Emanuel G. Kuflik, M.D., Adult & Child Dermatology Center,
Lakewood, NJ
Author of Cryosurgical Treatment for Skin Cancer*

"CRY-AC® is what dermatologists use every day in their practice.... I especially like Brymill's unit because it is very sturdy and holds up well after years of heavy use. They stand behind their units. Brymill is a company you can trust."

*Gloria F. Graham, M.D., Grahams' Dermatology Services, P.A.,
Winston-Salem, NC Author of cryosurgery chapter in R.
Wheeland's Cutaneous Surgery*

What other Physicians Say...

"Our family practice group has been using Brymill's CRY-AC® for about 15 years. We have achieved superb results with our patients. Brymill has constantly provided us with excellent service and a very personalized touch. I recommend the Brymill system without hesitation."

Neil H. Brooks, M.D., Past President, AAFP, Health Wise Medical Associates, Vernon, CT

"I use my Brymill CRY-AC® 20-30 times per week for treating everything from warts to molluscum to AK's to skin tags. It's incredibly easy to use. I trained my nurse practitioner to use it, and she uses it as often as I do."

Mark Moores, M.D., Webster Family Health Center, Dudley, MA

"Brymill's CRY-AC® is a very physician-friendly device for a busy family practitioner's office. All five doctors in our office use it, and at least one of the five of us uses it daily. The CRY-AC® quickly pays for itself many times over."

H. Andrew Selinger, M.D., Pro Health Physicians Family Medical Group, Bristol, CT